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User Manual

BodyTel™

The Telehealth Company

Soehnle Professional GmbH & Co. KG
Manfred-von-Ardenne-Allee 12
71522 Backnang
Germany

Vertrieb/Sales
BodyTel Europe GmbH
Hufelandstrasse 14
34537 Bad Wildungen
Germany

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Copyright ©September 2011 BodyTel Europe GmbH, all rights reserved.
No part of this manual may be reproduced or published in any form without the prior
written consent of BodyTel Europe GmbH.
BodyTel™, GlucoTel™, WeightTel™ and PressureTel™ are registered trademarks of
BodyTel Europe GmbH.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by BodyTel Europe GmbH under license. Other trademarks and trade names are those of their respective owners.

Java™ and all Java-based marks are trademarks or registered trademarks of Sun Microsystems, Inc. Other product and company names mentioned herein may be trademarks or trade names of their respective owners.

Android™ is a trademark of Google Inc. All rights reserved.

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
Table of Contents

1.	Functional elements	21
2.	General	22
2.1	Foreword	22
2.2	Scope of delivery	22
2.3	CE-mark	22
2.4	Manufacturer's declaration	22
3.	Introduction	24
3.1	General information on body weight	24
3.2	Purpose	24
4.	Using WeightTel	25
4.1	Before you start	25
4.2	Insert the batteries into WeightTel	25
4.3	Set up your BodyTel user account	26
4.4	Installing BodyTel Mobile	27
4.5	Pairing WeightTel	28
4.6	Positioning the WeightTel scale	29
4.7	Weighing operation with WeightTel	29
4.8	Transmitting data	30
4.9	Key to symbols	31
5.	Important instructions	32
5.1.	Safety measures	32
6.	Troubleshooting	33
6.1	Error codes	33
7.	Care and maintenance	34
7.1	Care	34
7.2	Measuring Technical Note	34
8.	Warranty and repair	35
8.1	Warranty	35
8.2	Repair	35
9.	Technical specifications	36
10.	EMV guidelines and Declaration of Conformity by Manufacturer	37

1 Functional elements

1. LCD-Display with 25 mm letters
2. Large 405 x 325 mm nonslip platform
3. Bluetooth® interface
4. Battery compartment
The battery compartment is located at the side of the scale
(see Chapter 4.2)
5. Bluetooth® pairing key
6. Bluetooth® status LED
7. Measurement unit switch (st,kg,lb)

2.1 Foreword

The integrated auto feedback system allows you to perform a first-class and comfortable weight measurement.

Measured values can be transmitted automatically via compatible mobile phones (smartphones) or the Home Gateway to your personal online patient file at the BodyTel Center (www.bodytel.com).

Thank you for your confidence. Should you have any service or product questions please contact us via E-Mail at support@bodytel.com or call the free BodyTel support at +800-BODYTEL (+800-0-2639835), Mo - Fri, 9 am - 6 pm.

2.2 Scope of delivery

1. WeightTel scale
2. User manual
3. Four AA batteries
4. Warranty card

2.3 CE-mark

WeightTel complies with the provisions of council directives 2004/108/EC and 2006/95/EC. As proof of conformity the scale carries the CE mark.

2.4 Manufacturer's declaration

Manufacturer: Soehnle Professional GmbH & Co. KG
Manfred-von-Ardenne-Allee 12
D-71522 Backnang

Backnang, Mai 2012, Soehnle Professional GmbH & Co. KG

2 General

Note

BodyTel Europe GmbH markets MEDISCALE BLUETOOTH under the name WeightTel.

Measuring accuracy

WeightTel measures with an accuracy of 0,35 % + 100g (1d)

Compatibility and costs

WeightTel can only be used with compatible mobile phones or the Home Gateway. A list of compatible devices can be found on www.bodytel.com. For transmitting measured values to your patient file an internet connection via LAN, WLAN or mobile internet is required. Please refer to your mobile phone plan for related costs.

3.1 General information on body weight

- The term “body weight” is used to describe the mass of a human body and is measured in kilograms throughout most countries of the world.
- There is no consent among health professionals as to what a person’s ideal weight should be.
- The World Health Organization (WHO) considers a person with a body mass index (BMI) below 18.5 as underweight. Underweight may lead to an insufficient supply of important nutrients such as proteins, fats, carbohydrates, vitamins and minerals. Severely underweight individuals may have poor stamina and a weak immune system, making them prone to infectious diseases.
- The World Health Organization defines obesity as follows: pre-obese with a BMI of 25 to 29.9, obesity grade I with a BMI of 30 to 34.9, obesity grade II with a BMI of 35 to 39.9, obesity class III with a BMI of 40 or higher. A raised BMI is linked to an increased risk of cardiovascular diseases, diabetes, musculoskeletal disorders and some cancers. Regular physical activity and healthier food choices (e. g. a limited intake of fats and carbohydrates) can help reduce obesity.
- Regular weight measurement can help you control your weight and achieve a healthier lifestyle.
- If you are in medical care because of your weight share your online patient file with your physician or take a print-out of the file with you. This will give your physician a detailed overview of your weight changes to support your therapy.

3.2 Purpose

- WeightTel is designed to measure body weight.
- Never change your medication on the basis of the obtained results. The evaluation of the measurements and the resulting treatment should only be carried out by your physician.
- Body weight measures should be taken without clothes or heavy accessories like mobile phones and wallets. Use similar conditions when taking multiple measurements to ensure comparable results.

4 Using WeightTel

4.1 Before you start


Before taking first measurements with WeightTel please insert the batteries into the scale (see chapter 4.2). Should you have received a preconfigured device (batteries included) please continue with chapter 4.3.

Note: To immediately transmit measured values to your online patient file activate your mobile phone's (smartphone) or Home Gateway's Bluetooth® interface and start the BodyTel Mobile application. Please keep the phone close to the WeightTel within a maximum range of 10 meters.

***Should you require assistance please contact the free support hotline:
+800-BODYTEL (+800-0-2639835), Mo - Fry, 9 am - 6 pm.***

4.2 Insert the batteries into WeightTel

Open the battery compartment located at the right side of the scale (seen from the back). Insert the four batteries that were supplied with the device. Please refer to the following picture (Abb. 2) for the correct orientation of the batteries.

4.3 Set up your BodyTel user account

To benefit from all the features of your new BodyTel telemedicine device please register for an user account at www.bodytel.com. This will enable you to automatically transmit weight measures from your WeightTel scale to your personal online patient file. Setting up is simple and free of charge. In order to register you will need:

- 1. A computer with an active internet connection**
- 2. Your e-mail address**

Step by Step instructions:

1. Open your internet browser and enter the following URL into the address field:
<http://www.bodytel.com>
2. Click on the "Create account" link, located in the upper right corner of the page.
3. In the following window enter your e-mail address into the corresponding field.
Click on "Register".
4. You will receive an e-mail confirmation to the specified address. This may take several minutes.
5. Open the e-mail sent from the address "server@bodytel.com" with the subject line "Email-Bestätigungs-Anfrage von www.bodytel.com".

Attention! Please check whether your e-mail client has moved the confirmation e-mail to the spam or junk mail folder. Should you have not received a confirmation e-mail within 30 Minutes after registering please repeat steps 1 through 5. Please contact our free service hotline at +800-BODYTEL (+800-0-2639835), Mo - Fry, 9 am - 6 pm, if the problem persists.

6. Click on the link beginning with "<https://secure.bodytel.com/de/optinconfirm/>".
7. A page will open in your browser reading "Welcome to the first visit to your BodyTel center". Please read and accept the "Terms of Service" and the "Privacy Policy" by checking the appropriate boxes. Click on "Next step".
8. Please fill out the required fields (sex, first name, last name) and choose a password. Click on "Next step".
9. Fill out the remaining fields and click on "Next step".
10. Choose your desired measurement units. The default settings should be: metric system, mg/dl, BE, 24 h, European.
11. Step 5 shows a summary of the data you have entered. Please check whether all details are correct and complete the registration procedure by clicking on "Done".
12. Your data has been saved and you will reach the personal front page "My BodyTel". This page is also displayed after you login at <http://www.bodytel.com>.

4 Using WeightTel

4.4 Installing BodyTel Mobile

Activate Bluetooth® on your mobile phone (smartphone) or Home Gateway and make sure visibility is set to "visible to all". There are two different ways to download the BodyTel Mobile application:

Option 1:

1. Open the internet browser on your mobile phone or Home Gateway and enter the following URL: <http://mobile.bodytel.com>.
2. On the BodyTel download page please choose your language by clicking on the respective link.
3. Please read and accept BodyTel's "Privacy Statement" and the "Business Terms" by checking the boxes. Click on "Send".
4. Please choose the version of the software compatible with your device.
5. The BodyTel Mobile software will now be downloaded. Please follow the instructions on your mobile phone or Home Gateway to install the software.

Option 2:

Download the software from Google's official Google Play Store. Please follow the instructions on your mobile phone or Home Gateway to install the software.

For further information please visit the BodyTel knowledge base at:

<http://www.bodytel-knowledgebase.com/DE>. (Only german language)

4 Using WeightTel

4.5 Pairing WeightTel

1. Please start the BodyTel Mobile application on your mobile phone (smartphone) or Home Gateway.
2. Enter your BodyTel user account data.
3. When starting the software for the first time the application will ask you whether you would like to pair a device now. Please confirm with "Yes".
4. Choose WeightTel from the device list.

Note: Your mobile phone or Home Gateway should be in close proximity to the WeightTel scale because the scale will always attempt to pair with the device sending the strongest Bluetooth® signal. Make sure that no other Bluetooth® enabled devices (mobile phones, computers, etc.) are close to the WeightTel while pairing.

5. Press the pairing key located at the bottom side of the scale's display (see left picture) until the display shows the word "COP". At the same time the Bluetooth® status LED will start flashing. The scale will now try to pair with the device sending the strongest Bluetooth® signal. The process is complete when the status LED is constantly on.

Note: If you press the pairing key during the pairing process the process will be terminated. To start over please repeat step 5.


6. Check your mobile phone or Home Gateway and confirm all connection requests.
7. Enter the Bluetooth® password "9676". On successful connection the status LED on the scale will flash twice and the display will show "PAI".

4 Using WeightTel

4.6 Positioning the WeightTel scale

Please note the following instructions regarding the correct positioning of the WeightTel scale before you proceed with the first measurements:


- The scale has to be placed on a flat, firm (uncarpeted) and non-slippery surface (check Fig. 3).
- The scale cannot perform correct measurements if it is not completely level.

4.7 Weighing operation with WeightTel

Note: Please read the introductory remarks concerning body weight and reduction in Chapter 3 "Introduction".

Place both feet firmly on the weighing platform. Hold your position until a weight is shown in the display. If you leave the platform before the measurement is completed WeightTel will automatically shut down.

4 Using WeightTel

If you want to transmit measured data into your online patient file at <http://www.bodytel.de> please remain on the scale for approximately 4 seconds after the weight has been displayed. Please follow the instructions given in Chapter 4.8 "Transmitting data". Should you not wish to transmit data leave the weighing platform. The scale will automatically shut down after approximately 5 seconds.

4.8 Transmitting data

During the transmission the measurements will be transmitted to your online patient file at <http://www.bodytel.com>.

During transmission the blue Bluetooth® status LED will start flashing. Upon completion, the LED will be activated constantly for 5 seconds. After that the scale will shut down. In case of a transmission failure the LED will turn off and the scale will display an error code (see chapter 6 "Troubleshooting") for 5 seconds.

Note: To transmit data the BodyTel Mobile software must be open and Bluetooth® must be activated on your mobile phone (smartphone) or Home Gateway.

4 Using WeightTel

4.9

Display

Key to symbols

LED

Operation/Situation

Bluetooth® pairing:

COP	rapid flashes	Bluetooth® pairing attempt
PAI	on	Strongest Bluetooth® device found
PAI	2 rapid flashes	Bluetooth® pairing successful

Weight Measurement:

Digits without unit (kg/lb)	off	Weight measurement initiated
Digits + unit (kg/lb)	off	Weight measurement complete

Bluetooth® data transmission:

Digits + unit (kg/lb)	slow flashes	Bluetooth® data transmission initiated
Digits + unit (kg/lb)	on	Bluetooth® data transmission successful
Alternating "rEtrY" and "1"	slow flashes	1st retry of Bluetooth® data transmission
Alternating "rEtrY" und "2"	slow flashes	2nd retry of Bluetooth® data transmission
noPAI	off	Data transmission failure, no pairing device found

Errors:

"Err" + error code	off	See chapter 6 "Troubleshooting"
lo	off	Replace batteries (see chapter 4.2)
overload	off	The maximum load (220 kg) is exceeded.

5 Important instructions



5.1 Safety measures

Attention:

Get the support of another person when using the scale if you feel frail or weak. Please note that a correct measurement can only be ensured if you stand on the scale unsupported during the actual weighing operation. Stop the procedure immediately whenever you feel dizzy or insecure.

In addition, please note:

- The scale is not water-resistant. Do not use the device in the bath tub or shower. If moisture has transgressed into the body of the scale remove the batteries. Stop using the scale and contact customer support.
- The maximum load is 220 kg.
- Consult your physician before you attempt to change your weight.
- Do not use the scale on slippery or wet surfaces.
- The scale is designed for indoor use only.

6 Troubleshooting

6.1 Error Codes

Error code	Description	Solution
Err 3	Internal communication error in Bluetooth® module, scale switches off after 5 seconds.	Please contact the vendor where you bought your BodyTel device.
Err 4	Bluetooth® transmission error, scale switches off after 5 seconds.	Please check whether a Bluetooth® device is in reach. Start the BodyTel Mobile app on your mobile phone (smartphone) or Home Gateway and check whether Bluetooth® connection is enabled.
Err 5	Bluetooth® transmission aborted after two attempts, scale switches off after 5 seconds.	Please check whether a Bluetooth® device is in reach. Start the BodyTel Mobile app on your mobile phone or Home Gateway and check whether Bluetooth® connection is enabled.
Err 8	Bluetooth® pairing refused, scale switches off after 5 seconds.	Please check whether you have entered the correct Bluetooth® password (9676) on your mobile phone or Home Gateway. Your mobile phone may be incompatible with WeightTel.

7 Care and maintenance

7.1 Care

Use a damp cloth and regular cleaning agents to clean the scale. Do not use any scouring agents.

7.2 Measuring Technical Note

This device is interference-free in accordance with the applicable EC Directive 2004/108/EC. The device is a measuring instrument. Wind, vibration, rapid temperature changes and sun exposure can have an effect on the measurement result.

Disposal of used electrical and electronic appliances including batteries:

According to our present knowledge the scale contains no environmentally harmful substances. This product is not to be treated as regular household waste, but should be handed in to an electrical/ electronic equipment recycling center. You can obtain further details from your local council, your municipal waste disposal company or the company from which you purchased the product.

Battery disposal: Do not dispose of batteries in the household trash. You will find these signs on batteries containing hazardous substances: Pb = battery contains lead, Cd = battery contains cadmium, Hg = battery contains mercury.

8 Warranty and repair

8.1 Warranty

1. BodyTel Europe GmbH grants a warranty of two years from date of sale. Proof of purchase must be supplied in the form of a correctly completed warranty card or a receipt.
2. Defects due to faulty material or workmanship will be repaired free of charge during the warranty period.
3. The warranty period is not extended as a result of warranty services being provided, either for the device as a whole or for any components.
4. The warranty excludes
 - a. any damage caused by improper use, e.g. by failure to observe the operating instructions.
 - b. damage caused by maintenance or interventions carried out by the purchaser or unauthorized third parties.
 - c. damage resulting from transport, either from the manufacturer to the customer or on the way to customer services.
 - d. component parts subject to normal wear and tear (e. g. batteries).
5. Liability is excluded for direct or indirect damage caused by the device even if the damage to the device falls within the terms of the warranty.
6. Unauthorized modification of the warranty card invalidates it.
7. Further claims for whatever reason are excluded.

BodyTel Europe GmbH does not warranty the batteries supplied.

8.2 Repair

In the case of a warranted repair please complete the warranty card giving precise details of the sender and the defect. Please contact our service department (see front cover) for shipping details.

9 Technical specifications

9.1 Technical specifications

Measuring range	max. 220 kg
Operating temperature	+5° C to +35° C (+41° F to 95° F)
Storage temperature	-20° C to +50° C (-4° F to +122° F)
Indicating range	max. 220 kg, 100 g steps
Accuracy	d = 100 g
Display elements	Display (LCD) Status LED
Operating elements	Bluetooth® pairing key Measurement unit switch
Data transmission	Bluetooth® Class I
Power supply	4 alkaline 1,5 V batteries (AA,LR6)
Unit dimensions (W x D x H)	325 mm x 480 mm x 120 mm
Tread surface	325 mm x 405 mm
Unit weight	4,9 kg

10 EMV guidelines and Declaration of Conformity by Manufacturer

10.1 EMV guidelines and Declaration of Conformity

The EMV guidelines and declaration of conformity for the MEDI SCALE BLUETOOTH scale with the product number 7841 (BodyTel product: WeightTel) are provided by the manufacturer on request or is available for download on the BodyTel website.

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Copyright ©September 2011 BodyTel Europe GmbH, alle Rechte vorbehalten.